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# COVID-19 Antigen Rapid Test Clinical Sensitivity and Specificity Study Report

## 1. Objective

The CLUNGENE<sup>®</sup> COVID-19 Antigen Rapid Test (hereinafter referred to as the CLUNGENE Device) manufactured by Hangzhou Clongene Biotech Co., Ltd. is a lateral flow immunoassay intended for the qualitative detection SARS-CoV-2 nucleocapsid antigens in nasopharyngeal swab and oropharyngeal swab from individuals who are suspected of COVID-19 by their healthcare provider.

This study is intended to evaluate the clinical performance, between the CLUNGENE Device and the comparator RT-PCR assay.

## 2. Method

A study of 770 direct nasopharyngeal swabs was performed. The specimens were prospectively collected from patients in community meeting Department of Public Health definition of a suspected COVID-19 case and being tested for SARS-CoV-2 part of routine medical care at 5 locations and tested at a single central laboratory.

Two nasopharyngeal swabs were collected from individual symptomatic (within 7 days of onset) or asymptomatic patients who were suspected of COVID-19. At all locations, one nasopharyngeal swab was tested directly with the COVID-19 Antigen Rapid Test according to product instructions for use, and the other swab was eluted in 3 mL viral transport media and tested with RT-PCR assay for detection of SARS-CoV-2. Swabs were randomly assigned to testing with the COVID-19 Antigen Rapid Test or RT-PCR assay and were tested by operators who were blinded to the RT-PCR test result.

The positive percent agreement (PPA) was calculated as  $100\% \times (\text{True Positive} / [\text{True Positive} + \text{False Negative}])$ . The negative percent agreement (NPA) was calculated as  $100\% \times (\text{True Negative} / [\text{True Negative} + \text{False Positive}])$ . Accuracy was calculated as  $100\% \times ([\text{True Positive} + \text{True Negative}] / \text{Total sample Qty})$ . The 95% (two-sided) confidence interval (CI) was calculated using the Wilson Score Method.

## 3. Comparator method

Real-Time Fluorescent RT-PCR Kit for Detecting SARS-CoV-2, manufactured by BGI Genomics Co. Ltd., is a real-time reverse transcription polymerase chain reaction (rRT-PCR) test. This product has got CE, NMPA certification and FDA Emergency Use Authorized. A specimen is positive for SARS-CoV-2 if the Ct value of ORF1ab gene is not higher than 37 and the Ct value of human housekeeping gene  $\beta$ -Actin is not higher than 35.

## 4. Enrollment criteria (inclusion/exclusion criteria)

### 4.1 Inclusion criteria

- Patients in community meeting Department of Public Health definition of a suspected COVID-19 case and being tested for SARS-CoV-2 part of routine medical care.
- Symptomatic (within 7 days of onset) or asymptomatic patients who were suspected of COVID-19.

### 4.2 Exclusion criteria

- Unable to obtain samples of information needed for the experiment
- Samples that have been contaminated or contaminated during sample storage
- Samples with inappropriate storage conditions

## 5. Result

The results are summarized in the following table.

The RT-PCR cycle threshold (Ct) is the relevant signal value. Lower Ct value indicate higher viral load. The

sensitivity was calculated for the different Ct value range (Ct value $\leq$ 33 and Ct value $\leq$ 37).

COVID-19Antigen		RT-PCR (Ct value $\leq$ 33)		Total
		Positive	Negative	
<b>CLUNGENE<sup>®</sup></b>	Positive	145	2	147
	Negative	3	593	596
Total		148	595	743

PPA (Ct $\leq$ 33): 98.0% (145/148), (95%CI: 94.2% ~99.3%)

NPA: 99.7% (593/595), (95%CI: 98.8% ~99.9%)

Accuracy: 99.3% ((145+593)/743), (95%CI: 98.4% ~99.7%)

COVID-19Antigen		RT-PCR (Ct value $\leq$ 37)		Total
		Positive	Negative	
<b>CLUNGENE<sup>®</sup></b>	Positive	161	2	163
	Negative	14	593	607
Total		175	595	770

PPA (Ct $\leq$ 37): 92.0% (161/175), (95%CI: 87.0% ~95.2%)

NPA: 99.7% (593/595), (95%CI: 98.8% ~99.9%)

Accuracy: 97.9% ((161+593)/770), (95%CI: 96.6% ~98.7%)

PPA - Positive Percent Agreement (Sensitivity)

NPA - Negative Percent Agreement (Specificity)

## 6. Conclusion

Taken together, the CLUNGENE<sup>®</sup> COVID-19 Antigen Rapid Test had a positive percent agreement (sensitivity) of 98.0% (95%CI: 94.2% ~99.3%) with specimens of a Ct count  $\leq$ 33, 92.0% (95%CI: 87.0% ~95.2%) with specimens of a Ct count  $\leq$ 37, negative percent agreement (specificity) of 99.7% (95%CI: 98.8% ~99.9%), and accuracy of 99.3% (95%CI: 98.4% ~99.7%) with specimens of a Ct count  $\leq$ 33, 97.9% (95%CI: 96.6% ~98.7%) with specimens of a Ct count  $\leq$ 37.

# COVID-19 Antigen Rapid Test Cassette (Nasal Swab) Clinical Sensitivity and Specificity Study Report

## 1. Objective

The CLUNGENE® COVID-19 Antigen Rapid Test Cassette (Nasal Swab) (hereinafter referred to as the CLUNGENE Device) manufactured by Hangzhou Clongene Biotech Co., Ltd. is a lateral flow immunoassay intended for the qualitative detection SARS-CoV-2 nucleocapsid antigens in nasal swab from individuals who are suspected of COVID-19 by their healthcare provider.

This study is intended to evaluate the clinical performance, between the CLUNGENE Device and the comparator RT-PCR assay.

## 2. Method

A study of 617 direct nasal swabs was performed. The specimens were prospectively collected from patients in community meeting Department of Public Health definition of a suspected COVID-19 case and being tested for SARS-CoV-2 part of routine medical care at 5 locations and tested at a single central laboratory.

Two nasal swabs were collected from individual symptomatic (within 7 days of onset) or asymptomatic patients who were suspected of COVID-19. At all locations, one nasal swab was tested directly with the COVID-19 Antigen Rapid Test Cassette (Nasal Swab) according to product instructions for use, and the other swab was eluted in 3 mL viral transport media and tested with RT-PCR assay for detection of SARS-CoV-2. Swabs were randomly assigned to testing with the COVID-19 Antigen Rapid Test Cassette (Nasal Swab) or RT-PCR assay and were tested by operators who were blinded to the RT-PCR test result.

The positive percent agreement (PPA) was calculated as  $100\% \times (\text{True Positive} / [\text{True Positive} + \text{False Negative}])$ . The negative percent agreement (NPA) was calculated as  $100\% \times (\text{True Negative} / [\text{True Negative} + \text{False Positive}])$ . Accuracy was calculated as  $100\% \times ([\text{True Positive} + \text{True Negative}] / \text{Total sample Qty})$ . The 95% (two-sided) confidence interval (CI) was calculated using the Wilson Score Method.

## 3. Comparator method

Real-Time Fluorescent RT-PCR Kit for Detecting SARS-CoV-2, manufactured by BGI Genomics Co. Ltd., is a real-time reverse transcription polymerase chain reaction (rRT-PCR) test. This product has got CE, NMPA certification and FDA Emergency Use Authorized. A specimen is positive for SARS-CoV-2 if the Ct value of ORF1ab gene is not higher than 37 and the Ct value of human housekeeping gene  $\beta$ -Actin is not higher than 35.

## 4. Enrollment criteria (inclusion/exclusion criteria)

### 4.1 Inclusion criteria

- Patients in community meeting Department of Public Health definition of a suspected COVID-19 case and being tested for SARS-CoV-2 part of routine medical care.
- Symptomatic (within 7 days of onset) or asymptomatic patients who were suspected of COVID-19.

### 4.2 Exclusion criteria

- Unable to obtain samples of information needed for the experiment
- Samples that have been contaminated or contaminated during sample storage
- Samples with inappropriate storage conditions

## 5. Result

The results are summarized in the following table.

The RT-PCR cycle threshold (Ct) is the relevant signal value. Lower Ct value indicate higher viral load. The sensitivity was calculated for the different Ct value range (Ct value  $\leq 30$  and Ct value  $\leq 37$ ).

COVID-19 Antigen		RT-PCR (Ct value $\leq$ 30)		Total
		Positive	Negative	
<b>CLUNGENE<sup>®</sup></b>	Positive	117	3	120
	Negative	3	462	465
Total		120	465	585

PPA (Ct $\leq$ 30):97.5% (117/120), (95%CI: 92.9%~99.2%)

NPA: 99.4% (462/465), (95%CI: 98.1%~99.8%)

COVID-19 Antigen		RT-PCR (Ct value $\leq$ 37)		Total
		Positive	Negative	
<b>CLUNGENE<sup>®</sup></b>	Positive	139	3	142
	Negative	13	462	475
Total		152	465	617

PPA (Ct $\leq$ 37):91.4% (139/152), (95%CI: 85.9%~94.9%)

NPA: 99.4% (462/465), (95%CI: 98.1%~99.8%)

PPA - Positive Percent Agreement (Sensitivity)

NPA - Negative Percent Agreement (Specificity)

## 6. Conclusion

Taken together, the CLUNGENE<sup>®</sup> COVID-19 Antigen Rapid Test Cassette (Nasal Swab) had a positive percent agreement (sensitivity) of 97.5% (95%CI: 92.9%~99.2%) with specimens of a Ct count  $\leq$ 30, 91.4% (95%CI: 85.9%~94.9%) with specimens of a Ct count  $\leq$ 37, negative percent agreement (specificity) of 99.4% (95%CI: 98.1%~99.8%).

## EC DECLARATION OF CONFORMITY

Name and address of the manufacturer: **Hangzhou Clongene Biotech Co., Ltd.  
No.1 Yichuang Road, Yuhang Sub-district  
Yuhang District  
311121 Hangzhou  
China**

We declare under our sole responsibility that

the medical device: **COVID-19 Antigen Rapid Test**

of class: **Other**  
according to article 9 of directive 98/79/EC

meets the provisions of the directive 98/79/EC and its transpositions in national laws which apply to it. The declaration is valid in connection with the "final inspection report" of the device.

Conformity assessment procedure: **Directive 98/79/EC Annex III**

Applicable standards: **EN ISO 13485:2016      EN ISO 15223-1:2016  
EN ISO 23640:2015      EN13612:2002/AC:2002  
EN 13975:2003      EN ISO 14971:2012  
EN ISO 18113-1:2011      EN ISO 18113-2:2011  
EN 62366-1:2015**

Name and address of the authorized representative: **Shanghai International Holding Corporation GmbH (Europe)  
Eiffestrasse 80  
20537 Hamburg  
Germany**

Hangzhou, July.15.2020

Place, date

Shujian Zheng, Legal representative

Name and function

# COVID-19 Antigen Rapid Test

The COVID-19 Antigen Rapid Test is a lateral flow immunoassay intended for the qualitative detection SARS-CoV-2 nucleocapsid antigens in nasopharyngeal swab, nasal swab or oropharyngeal swab from individuals who are suspected of COVID-19 by their healthcare provider.

## Features

CE

Easy to collect samples.

No equipment required.

Instant result at 15 minutes.

Results are clearly visible.


Suitable for large-scale rapid screening.

## Components


Extraction Reagent  
Tube x 25

Sterilized Swab  
x 25

Test Cassette  
x 25

Work Station  
x 1

Package Insert  
x 1

## Samples

Nasopharyngeal Swab


Nasal Swab

Oropharyngeal Swab

## Test Procedure

## Performance Characteristics

### Clinical Performance

#### Nasopharyngeal Swab

770 nasopharyngeal swabs were collected from individual symptomatic patients by using CLUNGENE® COVID-19 Antigen Rapid Test. The RT-PCR cycle threshold (Ct) is the relevant signal value. Lower Ct value indicate higher viral load. The sensitivity was calculated for the different Ct value range (Ct value  $\leq 33$  and Ct value  $\leq 37$ ).

Summary data as below:

COVID-19 Antigen		RT-PCR (Ct value $\leq 33$ )		Total
		Positive	Negative	
CLUNGENE®	Positive	145	2	147
	Negative	3	593	596
Total		148	595	743

PPA (Ct  $\leq 33$ ): 98.0% (145/148), (95%CI: 94.2%~99.3%)  
 NPA: 99.7% (593/595), (95%CI: 98.8%~99.9%)

COVID-19 Antigen		RT-PCR (Ct value $\leq 37$ )		Total
		Positive	Negative	
CLUNGENE®	Positive	161	2	163
	Negative	14	593	607
Total		175	595	770

PPA (Ct  $\leq 37$ ): 92.0% (161/175), (95%CI: 87.0%~95.2%)  
 NPA: 99.7% (593/595), (95%CI: 98.8%~99.9%)

PPA - Positive Percent Agreement (Sensitivity)

NPA - Negative Percent Agreement (Specificity)

## Nasal Swab

617 nasal swabs were collected from individual symptomatic patients by using CLUNGENE® COVID-19 Antigen Rapid Test. The sensitivity was calculated for the different Ct value range (Ct value $\leq$ 33 and Ct value $\leq$ 37).

Summary data as below:

COVID-19 Antigen		RT-PCR (Ct value $\leq$ 33)		Total
		Positive	Negative	
CLUNGENE®	Positive	132	3	135
	Negative	4	462	466
Total		136	465	601

PPA (Ct  $\leq$  33): 97.1% (132/136), (95%CI: 92.7%~98.9%)

NPA: 99.4% (462/465), (95%CI: 98.1%~99.8%)

COVID-19 Antigen		RT-PCR (Ct value $\leq$ 37)		Total
		Positive	Negative	
CLUNGENE®	Positive	139	3	142
	Negative	13	462	475
Total		152	465	617

PPA (Ct  $\leq$  37): 91.4% (139/152), (95%CI: 85.9%~94.9%)

NPA: 99.4% (462/465), (95%CI: 98.1%~99.8%)

PPA - Positive Percent Agreement (Sensitivity)

NPA - Negative Percent Agreement (Specificity)

### Limit of Detection (Analytical Sensitivity)

The study used cultured SARS-CoV-2 virus which is heat inactivated and spiked into nasopharyngeal swab specimen. The Limit of Detection (LoD) is  $5.7 \times 10^2$  TCID<sub>50</sub>/mL.

### Cross Reactivity (Analytical Specificity)

32 commensal and pathogenic microorganisms that may be present in the nasal cavity have been evaluated and no cross-reactivity was observed.

### Interference

17 potential interference substances with different concentration were evaluated and found no affect to the test performance.

### High-dose Hook Effect

The COVID-19 Antigen Rapid Test was tested up to  $1.0 \times 10^{5.67}$  TCID<sub>50</sub>/mL of inactivated SARS-CoV-2 and no high-dose hook effect was observed.

**Hangzhou Clongene Biotech Co., Ltd.**

Add: No.1 Yichuang Road, Yuhang Sub-district, Yuhang District, 311121 Hangzhou, China

Tel: +86 (0)571-88262120

Fax: +86 (0)571-88261752

## EC DECLARATION OF CONFORMITY

Name and address of the manufacturer: **Hangzhou Clongene Biotech Co., Ltd.**  
**No.1 Yichuang Road, Yuhang Sub-district**  
**Yuhang District**  
**311121 Hangzhou**  
**China**

We declare under our sole responsibility that

the medical device: **COVID-19 Antigen Rapid Test Cassette (Nasal Swab)**

of class: **Other**  
according to article 9 of directive 98/79/EC

meets the provisions of the directive 98/79/EC and its transpositions in national laws which apply to it. The declaration is valid in connection with the "final inspection report" of the device.

Conformity assessment procedure: **Directive 98/79/EC Annex III, excluding Section 6**

Applicable standards: **EN ISO 13485:2016**                      **EN ISO 15223-1:2016**  
**EN ISO 18113-1:2011**                      **EN ISO 18113-2:2011**  
**EN ISO 23640:2015**                      **EN ISO 14971:2019**  
**EN 13641:2002**                              **EN13612:2002/AC:2002**  
**EN 13975:2003**                              **EN 62366-1:2015**  
**ISO 17511:2020**

Name and address of the authorized representative: **Shanghai International Holding Corporation GmbH (Europe)**  
**Eiffestrasse 80**  
**20537 Hamburg**  
**Germany**

Hangzhou, February 04, 2021.

Place, date

Shujian Zheng, Legal representative

Name and function

# Certificate

The Certification Body of  
TÜV Rheinland LGA Products GmbH

hereby certifies that the organization  
**Hangzhou Clongene Biotech  
Co., Ltd.**  
**No. 1 Yichuang Road, Yuhang Sub-district  
Yuhang District  
311121 Hangzhou  
P.R. China**

has established and applies a quality management system for medical devices  
for the following scope:

**Design and Development, Manufacture and Distribution of  
In-vitro Diagnostic Rapid Test of Fertility, Drug of Abuse,  
Infectious Diseases, Tumour Markers and Cardiac Markers**

Proof has been furnished that the requirements specified in

**EN ISO 13485:2016**

are fulfilled. The quality management system is subject to yearly surveillance.

Effective Date: 2020-11-25  
Certificate Registration No.: SX 60152722 0001  
An audit was performed. Report No.: 15073650 008  
This Certificate is valid until: 2023-11-18

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July 6, 2021

To whom it may concern,

We, as the manufacturer of the CLUNGENE COVID-19 Antigen Rapid Test, assure that the following major SARS-CoV-2 variants which seem to spread more easily and quickly than other variants and are circulating globally, can be detected by CLUNGENE COVID-19 Antigen Rapid Test:

- ✓ **B.1.1.7 (Alpha)**: This variant was initially detected in the UK;
- ✓ **B.1.351 (Beta)**: This variant was initially detected in South Africa in December 2020;
- ✓ **P.1 (Gamma)**: This variant was initially identified in travelers from Brazil, who were tested during routine screening at an airport in Japan, in early January;
- ✓ **B.1.617.2 (Delta)**: This variant was initially detected in India.

The recombinant N proteins of these variants can be detected by the CLUNGENE COVID-19 Antigen Rapid Test, at concentrations same as that of the ordinary recombinant N protein without variation (amino acid: Met1-Ala419, accession # QHO62884.1).

Yours sincerely,

Shujian Zheng,

Managing Director of Hangzhou Clongene Biotech Co., Ltd.

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EC REP Shanghai International Holding Corp.GmbH (Europe)