

An official EU website



[Live, work, travel in the EU](#)

## COVID-19 In Vitro Diagnostic Devices and Test Methods Database

[Home](#) > [COVID-19 In Vitro Diagnostic Medical Devices](#) >

**COVID-19 In Vitro Diagnostic Medical Device - detail**

---

# COVID-19 In Vitro Diagnostic Medical Device - detail

## SARS-CoV-2 Antigen Kit (Colloidal Gold)

**Manufactured by Goldsite Diagnostics Inc, China -**

[en.goldsite.com.cn/](https://en.goldsite.com.cn/) 

Device

identification

number

1197

CE Marking ✓ Yes

HSC ✓ Yes

common list

HSC mutual ✓ Yes

recognition

Format Near POC / POC

Target Antigen

Specimen Nasal swab, Other

Commercial Commercialised  
Status

Last Update 2021-07-13 06:05:32 CET

## Comments

[Show HSC list status history](#) ▾

---

Assay Type	Immuno-Antigen
Rapid Diagnostic	Yes
Reader Required	No
Method	Immunoassay
Measurement	
Qualitative	
Time	12 minutes
Detection Principle	Colloidal gold

---

The database contains publicly available In Vitro Diagnostic Medical Devices for COVID-19 and it is being updated periodically. Please note that additional performance (as retrieved from manufacturers web pages) is provided only for devices commercially available with CE-IVD mark.

[Acknowledgements](#)

## COVID-19 Test Methods and Devices

This site is managed by the Joint Research Centre

[Contact us](#)

[Acknowledgements](#)

---

More information on:

[Live, work, travel in the EU](#)

---

## European Commission

[Contact the European Commission](#)

[Follow the European Commission on social media](#)

[Resources for partners](#)

[Language policy](#)

[Cookies](#)

[Privacy policy](#)

[Legal notice](#)